



# Test Definition: FMPRE

Maple Red (Acer rubrum) IgE

## Overview

### Method Name

Enzyme Immunoassay (EIA)

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

0.5 mL

### Reject Due To

Hemolysis:	NA
Thawing:	Warm OK; Cold OK
Lipemia:	Mild OK; Gross reject
Icterus:	NA
Other:	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

## Clinical & Interpretive

### Reference Values

<0.35 kU/L

## Interpretation

<a href="#">Class</a>	<a href="#">IgE (kU/L)</a>	<a href="#">Comment</a>
0	<0.35	Below Detection
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>=50	Very Strong Positive

## Performance

### Method Description

The conventional EIA uses allergen-coated discs from several suppliers and an enzyme-labeled anti-IgE.

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

2 to 7 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

[This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.](#)

### CPT Code Information

86003

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FMPRE	Maple Red (Acer rubrum) IgE	7479-9

Result ID	Test Result Name	Result LOINC® Value
Z3094	Maple Red (Acer rubrum) IgE	7479-9
Z3095	Class	102427-2