



Test Definition: FFRFT

Anti-Synthetase Profile

Overview

Method Name

RIPA Gel Radiography; Enzyme-Linked Immunosorbent Assay (ELISA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Collection Container/Tube: Red top or serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL Serum

Collection Instructions:

1. Within 1 hour of collection, centrifuge and aliquot 3 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume

Serum: 1 mL; Note: This volume does not allow for repeat testing.

Reject Due To

Gross hemolysis:	Reject
Lipemia:	Reject
Icterus	Reject
Anything other than serum	Reject
Bacterial contamination	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	7 days	
	Frozen	60 days	

Clinical & Interpretive

Clinical Information

Refer to www.labcorp.com/test-menu

Reference Values

Anti-PL-7 Ab, Anti-PL-12 Ab, Anti-EJ Ab, Anti-OJ Ab: Negative

Anti-Jo-1 Ab:

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

16 to 25 days

Performing Laboratory Location

Esoterix Endocrinology

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

CPT Code Information

83516 x 4

86235

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFRFT	Anti-Synthetase Profile	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5782	Anti-Jo-1 Ab	35333-4
Z5783	Anti-PL-7 Ab	33772-5
Z5837	Anti-PL-12 Ab	33771-7
Z5823	Anti-EJ Ab	45149-2
Z5824	Anti-OJ Ab	45152-6