



# Test Definition: FADFL

Adenosine Deaminase in Peritoneal Fluid

## Overview

### Method Name

Quantitative Spectrophotometry

### NY State Available

Yes

## Specimen

### Specimen Type

Peritoneal

### Specimen Required

**Specimen Type:** Peritoneal fluid (ascites, abdominal fluid, paracentesis fluid)

**Collection Container/Tube:** Leak-proof container

**Submission Container/Tube:** Plastic vial

**Specimen volume:** 0.5 mL

#### Collection Instructions:

1. Centrifuge specimen at ambient temperature.
2. Transfer 0.5 mL into a plastic vial and freeze
3. Send frozen. **Specimen must remain frozen** until received at the performing lab.
4. Indicate source on requisition.

### Specimen Minimum Volume

0.2 mL

### Reject Due To

|                              |        |
|------------------------------|--------|
| Bronchoalveolar lavage (BAL) | Reject |
| Glass tubes                  | Reject |
| Turbid specimens             | Reject |
| Whole blood                  | Reject |

### Specimen Stability Information

| Specimen Type | Temperature        | Time    | Special Container |
|---------------|--------------------|---------|-------------------|
| Peritoneal    | Frozen (preferred) | 30 days |                   |
|               | Refrigerated       | 7 days  |                   |

## Clinical & Interpretive

### Clinical Information

Refer to <https://ltd.aruplab.com/>

### Reference Values

0-30 U/L

## Performance

### PDF Report

No

### Day(s) Performed

Tuesday, Thursday, Sunday

### Report Available

1 to 8 days

### Performing Laboratory Location

ARUP Laboratories

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

### CPT Code Information

84311

### LOINC® Information

| Test ID | Test Order Name                    | Order LOINC® Value |
|---------|------------------------------------|--------------------|
| FADFL   | Adenosine Deaminase Peritoneal Fld | 49759-4            |

# Test Definition: FADFL

Adenosine Deaminase in Peritoneal Fluid

| Result ID | Test Result Name                   | Result LOINC® Value |
|-----------|------------------------------------|---------------------|
| FADFL     | Adenosine Deaminase Peritoneal Fld | 49759-4             |