



Test Definition: FCSC

Synthetic Cannabinoid Metabolites Screen,
Expanded, Urine

Overview

Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FCSCU	Synthetic Cannabinoid Confirm, U	No	No

Testing Algorithm

Confirmation test, FCSCU, will be performed, if appropriate, at no additional charge.

Method Name

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Urine

Specimen Required

*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Specimen Type: Urine

Container/Tube: Plastic preservative-free urine container

Specimen Volume: 5 mL

Collection Instructions: Collect 5 mL random urine specimen without preservative. Send specimen refrigerated in a plastic preservative-free urine container.

Note: Known interference(s): FUBINACA 3,3-dimethylbutanoic acid: Quetiapine

Specimen Minimum Volume

2.4 mL

Reject Due To

Thawing:	Warm OK; Cold OK
----------	------------------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	30 days	
	Ambient	30 days	
	Frozen	180 days	

Clinical & Interpretive

Reference Values

Reporting limit determined each analysis.

<u>Test</u>	<u>Result</u>	<u>Units</u>
5F-PINACA 3-methylbutanoic acid	None Detected	ng/mL
4F-BINACA 3,3-dimethylbutanoic acid	None Detected	ng/mL
5F-PICA 3,3-dimethylbutanoic acid	None Detected	ng/mL
5F-PINACA 3,3-dimethylbutanoic acid	None Detected	ng/mL
MDMB-4en-PINACA butanoic acid	None Detected	ng/mL
FUBINACA 3-methylbutanoic acid	None Detected	ng/mL
FUBINACA 3,3-dimethylbutanoic acid	None Detected	ng/mL
4-carboxy-NA-PIM	None Detected	ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

10 to 16 days

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.

- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80307

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCSC	Synthetic Cannabinoid Screen, Urine	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6094	5F-PINACA 3-methylbutanoic acid	Not Provided
Z6095	4F-BINACA 3,3-dimethylbutanoic acid	Not Provided
Z6096	5F-PICA 3,3-dimethylbutanoic acid	Not Provided
Z6097	5F-PINACA 3,3-dimethylbutanoic acid	Not Provided
Z6098	MDMB-4en-PINACA butanoic acid	Not Provided
Z6099	FUBINACA 3-methylbutanoic acid	Not Provided
Z6100	FUBINACA 3,3-dimethylbutanoic acid	Not Provided
Z6101	4-carboxy-NA-PIM	Not Provided