



# Test Definition: THEO

Theophylline, Serum

## Overview

### Useful For

Assessing and adjusting theophylline dosage for optimal therapeutic level

Assessing theophylline toxicity

### Method Name

Kinetic Interaction of Microparticles in a Solution (KIMS)

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.5 mL

**Collection Instructions:**

1. Serum gel tubes should be centrifuged within 2 hours of collection.
2. Red-top tubes should be centrifuged and the serum aliquoted into a plastic vial within 2 hours of collection.

### Forms

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

### Specimen Minimum Volume

0.25 mL

### Reject Due To

Gross hemolysis	Reject
-----------------	--------

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	28 days	

## Clinical & Interpretive

### Clinical Information

Theophylline and its congener, aminophylline, are used to relax smooth muscles of the bronchial airways and pulmonary blood vessels to relieve and prevent symptoms of asthma and bronchospasm.

Theophylline is typically administered orally at a dose of 400 mg/day or 6 mg/kg, whichever is lower, or intravenously as aminophylline at 0.4 mg/kg/hour. Oral dosage may be increased at 200 mg increments to a maximum of 900 mg/day, or 13 mg/kg if the steady-state blood concentration is within the therapeutic range of 8.0 to 20.0 mcg/mL.

Theophylline has a half-life of approximately 4 hours in children and adult smokers, and 8.7 hours in nonsmoking adults. The volume of distribution is approximately 0.5 L/kg, and the drug is approximately 40% protein bound. Theophylline exhibits zero-order clearance kinetics like phenytoin, small increases in dose yield disproportionately large increases in blood concentration.

Coadministration of cimetidine and erythromycin will significantly inhibit theophylline clearance, requiring dosage reduction. Other drugs such as allopurinol, ciprofloxacin, oral contraceptives, and propranolol inhibit theophylline clearance to a lesser degree.

Smoking induces the synthesis of cytochrome P448, the antipyrine-dependent cytochrome, which significantly increases the rate of metabolism of theophylline. Drugs such as phenobarbital, phenytoin, carbamazepine, and rifampin slightly increase the rate at which the drug is cleared.

Theophylline exhibits rather severe toxicity that is proportional to blood level.

### Reference Values

Therapeutic:

Bronchodilation: 8.0-20.0 mcg/mL

Neonatal apnea (< or =4 weeks old): 6.0-13.0 mcg/mL

Critical value: >20.0 mcg/mL

### Interpretation

Response to theophylline is directly proportional to the serum level.

Patients usually receive the best response when the serum level is above 8.0 mcg/mL, with minimal toxicity experienced as long as the level is less than or equal to 20.0 mcg/mL.

### Cautions

Many drugs affect the plasma level of this drug (as outlined in Clinical Information).

---

**Clinical Reference**

1. Milone MC, Shaw LM. Therapeutic drugs and their management. In: Rifai N, Chiu RWK, Young I, Burnham C-AD, Wittwer CT, eds. Tietz Textbook of Laboratory Medicine. 7th ed. Elsevier; 2023:420-453.e9
2. Brunton LL, Hilal-Dandan R, Knollmann BC, eds. Goodman and Gilman's: The Pharmacological Basis of Therapeutics. McGraw-Hill; 2018

**Performance****Method Description**

The assay is based on the kinetic interaction of microparticles in a solution. Theophylline antibody is covalently coupled to microparticles and the drug derivative is linked to a macromolecule. The kinetic interaction of microparticles in solutions is induced by binding of drug-conjugate to the antibody on the microparticles and is inhibited by the presence of theophylline in the sample. A competitive reaction takes place between the drug conjugate and theophylline in the serum sample for binding to the theophylline antibody on the microparticles. The resulting kinetic interaction of microparticles is indirectly proportional to the amount of drug present in the sample. (Package insert: Theophylline reagent. Roche Diagnostics; 09/2017)

**PDF Report**

No

**Day(s) Performed**

Monday through Saturday

**Report Available**

Same day/1 to 2 days

**Specimen Retention Time**

2 weeks

**Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Superior Drive

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been modified from the manufacturer's instructions. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food

---

and Drug Administration.

**CPT Code Information**

80198

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
THEO	Theophylline, S	4049-3

Result ID	Test Result Name	Result LOINC® Value
8661	Theophylline, S	4049-3