



# Test Definition: TLCU

Immunoglobulin Total Light Chains, Urine

## Overview

### Useful For

Monitoring patients whose urine demonstrates large M-spikes

Confirming the quantitation of specimens that show M-spikes by electrophoresis

Detecting urine monoclonal proteins and identification of specimens that need urine protein electrophoresis

### Profile Information

Test Id	Reporting Name	Available Separately	Always Performed
KTLCU	Kappa Total Light Chain, U	No	Yes
LTLCU	Lambda Total Light Chain, U	No	Yes
KLTRU	Kappa/Lambda TLC Ratio, U	No	Yes

### Special Instructions

- [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#)

### Method Name

Nephelometry

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Ordering Guidance

If serum is being submitted on the same patient for FLCS / Immunoglobulin Free Light Chains, Serum; that test should be ordered under a different order number.

### Specimen Required

Submit only 1 of the following specimens:

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Specimen Type:** Random urine

**Collection Container/Tube:** Clean, plastic urine collection container

**Submission Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 1 mL

**Collection Instructions:**

1. Collect a random urine specimen.
2. Refrigerate after collection and send refrigerate.

**Specimen Type:** 24-Hour urine

**Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 1 mL

**Collection Instructions:**

1. Collect urine for 24 hours.
2. No additive needed for preservation.
3. Urine may be kept ambient during the collection period but should be refrigerated within 4 hours of collection completion.

**Additional Information:** See [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#) for multiple collections.

**Urine Preservative Collection Options**

**Note:** The application of temperature controls **must occur within 4 hours of completion** of the collection.

Ambient (No additive)	OK <72 hours
Refrigerate (No additive)	Preferred
Frozen (No additive)	OK
50% Acetic Acid	No
Boric Acid	No
Diazolidinyl Urea	No
6M Hydrochloric Acid	No
6M Nitric Acid	No
Sodium Carbonate	No
Thymol	No
Toluene	No

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	7 days	

	Ambient	72 hours	
	Frozen	20 days	

## Clinical & Interpretive

### Clinical Information

Immunoglobulin light chains are usually cleared from blood through the renal glomeruli and reabsorbed in the proximal tubules so that urine light-chain concentrations are very low or undetectable. The production of large amounts of monoclonal light chains, however, can overwhelm this reabsorption mechanism. The detection of monoclonal light chains in the urine (Bence Jones proteinuria) has been used as a diagnostic marker for multiple myeloma since the report by Dr. H. Bence Jones in 1847.

Current laboratory procedures employ protein electrophoresis and isotype testing for the identification and characterization of urine monoclonal light chains, which may be present in large enough amounts to also be quantitated as an M-spike on protein electrophoresis. The electrophoretic M-spike is the recommended method of monitoring monoclonal gammopathies, such as multiple myeloma. Monitoring the urine M-spike is especially useful in patients with light-chain multiple myeloma in whom the serum M-spike is very small or absent, but the urine M-spike is large.

Just as quantitative serum immunoglobulins by immunonephelometry are a complement to M-spike quantitation by serum electrophoresis, this quantitative urine light-chain assay may be used to complement urine M-spike quantitation by electrophoresis.

### Reference Values

KAPPA TOTAL LIGHT CHAIN

<0.9 mg/dL

LAMBDA TOTAL LIGHT CHAIN

<0.7 mg/dL

KAPPA/LAMBDA RATIO

0.7-6.2

### Interpretation

A kappa/lambda (K/L) ratio greater than 6.2 suggests the presence of monoclonal kappa light chains.

A K/L ratio less than 0.7 suggests the presence of monoclonal lambda light chains.

In 24-hour specimens, a greater than 90% increase in concentration suggests progression or relapse; a greater than 90% decrease suggests treatment response.

Increased kappa and/or lambda light chains may be seen in benign (polyclonal) and neoplastic (monoclonal) disorders.

### Cautions

Unlike the electrophoretic M-spike, this immunoassay quantitates both polyclonal and monoclonal light chains and is therefore not sensitive for detecting small monoclonal abnormalities. A normal kappa/lambda (K/L) ratio does not rule

out a monoclonal protein, and an abnormal ratio does not identify a monoclonal protein. Urine protein electrophoresis and isotype testing are more sensitive and specific.

The quantitation of urine kappa light chain by immunonephelometry yields results that are approximately 2 times the values from the electrophoresis M-spike. Sequential results should be compared to previous results obtained by the same methodology.

### Supportive Data

In a study of 168 urine samples with a monoclonal light chain detected by immunofixation electrophoresis (IFE), there were 20 samples with a normal kappa/lambda (K/L) ratio. These samples had either no M-spike (n=13) or M-spikes <0.5 mg/dL. Conversely, among the 148 cases with an abnormal K/L ratio, there were 12 samples with no M-spike indicating that there is no clear M-spike value at which the K/L ratio identifies monoclonal light chains. In patients with an M-spike, the relationship between the kappa and lambda light-chain quantitation and the size of the M-spike had good correlation (kappa,  $r[2]=0.94$ ; lambda,  $r[2]=0.71$ ) and the regression lines had slopes of 2.4 of kappa and 1.1 for lambda.

Interestingly, there was a single case in which the K/L ratio was 24 and the free light-chain K/L ratio was 58, but the IFE showed polyclonal light chains. The patient was post-transplant for a kappa light-chain multiple myeloma and presumably had multiple forms of a monoclonal kappa light chain that migrated in a smear and was a false-negative by IFE.

### Clinical Reference

1. Kumar S, Paiva B, Anderson KC, et al. International Myeloma Working Group consensus criteria for response and minimal residual disease assessment in multiple myeloma. *Lancet Oncol.* 2016;17(8):e328-e346
2. Leung N, Barnidge DR, Hutchison CA. Laboratory testing in monoclonal gammopathy of renal significance (MGRS). *Clin Chem Lab Med.* 2016;54(6):929-937
3. Willrich MA, Katzmann JA. Laboratory testing requirements for diagnosis and follow-up of multiple myeloma and related plasma cell dyscrasias. *Clin Chem Lab Med.* 2016;54(6):907-919
4. Rajkumar SV, Kyle RA. Multiple myeloma: diagnosis and treatment. *Mayo Clin Proc.* 2005;80(10):1371-1382

### Performance

#### Method Description

In this Siemens Nephelometer II method, the light scattered by the antigen-antibody complexes is measured. The intensity of the measured scattered light is proportional to the amount of antigen-antibody complexes in the sample under certain conditions. If the antibody volume is kept constant, the signal behaves proportionally to the antigen volume.

A reference curve is generated by a standard with a known antigen content on which the scattered light signals of the samples can be evaluated and calculated as an antigen concentration. Antigen-antibody complexes are formed when a sample containing antigen and the corresponding antiserum are put into a cuvette. A light beam is generated with a light emitting diode (LED), which is transmitted through the cuvette. The light is scattered onto the immuno-complexes that are present. Antigen and antibody are mixed in the initial measurement, but no complex is yet formed. An antigen-antibody complex is formed in the final measurement.

The result is calculated by subtracting value of the final measurement from the initial measurement. The distribution of intensity of the scattered light depends on the ratio of the particle size of the antigen-antibody complexes to the radiated wavelength. (Instruction manual: Siemens Nephelometer II. Siemens, Inc; Version 2.4, 07/2019; Addendum to the Instruction Manual 2.3, 08/2017)

**PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

Same day/1 to 3 days

**Specimen Retention Time**

14 days

**Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Superior Drive

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been modified from the manufacturer's instructions. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

83883 x 2

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
TLCU	Immunoglobulin Total Light Chains,U	44792-0

  

Result ID	Test Result Name	Result LOINC® Value
KLTRU	Kappa/Lambda TLC Ratio, U	33559-6
KTLCU	Kappa Total Light Chain, U	27365-6
LTLCU	Lambda Total Light Chain, U	27394-6