



# Test Definition: FFTLB

Toluene, Occupational Exposure, Blood

## Overview

### Method Name

Headspace Gas Chromatography/Flame Ionization Detection (GC-FID)

### NY State Available

Yes

## Specimen

### Specimen Type

WB Sodium Heparin

### Specimen Required

Draw blood in 2 green top (sodium heparin) tubes. Tubes should be filled to prevent loss of volatile compound into headspace. Send 20 mL of sodium heparin whole blood refrigerated.

### Specimen Minimum Volume

5 mL

### Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
WB Sodium Heparin	Refrigerated (preferred)	14 days	
	Frozen	365 days	

## Clinical & Interpretive

### Reference Values

Units: mg/L

Normal (unexposed population): None detected

**Exposed:**

Biological Exposure Index (BEI): 0.05 mg/L (prior to last shift of workweek)

Biological Tolerance Value (BAT): 1.0 mg/L (end of exposure or end of shift)

**Toxic:**

Blood levels between 50 and 79 mg/L were found in people who died of acute toluene inhalation.

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

7 to 11 days

**Specimen Retention Time**

2 weeks

**Performing Laboratory Location**

Medtox Laboratories, Inc.

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**CPT Code Information**

84600

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FFTLB	Toluene, Blood	9637-0

Result ID	Test Result Name	Result LOINC® Value
-----------	------------------	---------------------

---

Z1203	Toluene	9637-0
-------	---------	--------