



# Test Definition: FPHET

## Anti-Phosphatidylethanolamine Panel

### Overview

#### Method Name

Enzyme-Linked Immunosorbent Assay (ELISA)

#### NY State Available

Yes

### Specimen

#### Specimen Type

Serum

#### Specimen Required

##### Collection Container/Tube:

**Preferred:** Red top

**Acceptable:** Serum gel

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 3 mL Serum

##### Collection Instructions:

1. Allow blood to clot for 30 minutes.
1. Centrifuge and aliquot 3 mL of serum into a plastic vial.
2. Send frozen.

#### Specimen Minimum Volume

Serum: 0.5 mL

#### Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Gross icterus	Reject
Specimens containing microbial contamination	Reject

#### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Serum	Frozen (preferred)	30 days	
	Refrigerated	14 days	

## Clinical & Interpretive

### Clinical Information

Refer to [www.bioagilytix.com/solutions/capabilities/diagnostic-testing/](http://www.bioagilytix.com/solutions/capabilities/diagnostic-testing/)

### Reference Values

Anti-Phosphatidylethanolamine IgA  
<12.0 U/mL

Anti-Phosphatidylethanolamine IgG  
<12.0 U/mL

Anti-Phosphatidylethanolamine IgM  
<12.0 u/mL

Reference range applies to anti-phosphatidylethanolamine IgA, IgG, & IgM

Normal: <12.0

Equivocal: 12.0-18.0

Elevated: >18.0

## Performance

### PDF Report

No

### Day(s) Performed

Wednesday

### Report Available

3 to 11 days

### Performing Laboratory Location

BioAgilytix Diagnostics

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

The performance characteristics of the listed assays were validated by BioAgilytix Diagnostics. The US FDA has not approved or cleared these tests. The results of these assays can be used for clinical diagnosis without FDA approval. BioAgilytix Diagnostics is a CLIA certified, CAP accredited laboratory for performing high complexity assays.

### CPT Code Information

83520 x 3

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPHET	Anti-Phosphatidylethanolamine Panel	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z0150	Anti-Phosphatidylethanolamine IgG	13076-5
Z0142	Anti-Phosphatidylethanolamine IgM	13077-3
Z0143	Anti-Phosphatidylethanolamine IgA	13078-1