

Instructions: This document must be completed with blue or black ink. If any other colored ink is used, this Chain of Custody process will be rejected.

Client Information (required)

Client Name		
Client Account No.		
Client Phone		
Street Address		
City	State	ZIP Code

Patient Information (required)

Patient ID (Medical Record No.)	
Patient Name (Last, First Middle)	
Sex <input type="checkbox"/> Male <input type="checkbox"/> Female	Birth Date (mm-dd-yyyy)
Collection Date (mm-dd-yyyy)	Time <input type="checkbox"/> am <input type="checkbox"/> pm

Submitting Healthcare Professional Information (required)

Submitting/Referring Healthcare Professional Name (Last, First)

Fill in only if Call Back is required.

Phone (with area code)	Fax* (with area code)
National Provider Identification (NPI)	

**Fax number given must be from a fax machine that complies with applicable HIPAA regulation.*

MCL Internal Use Only

Both pages of the original form must accompany the specimen to Mayo Clinic Laboratories. Photocopies and faxes are not valid for Chain of Custody testing. It is the client's responsibility to maintain documentation of the order.

All items in this section below must be completed or Chain of Custody will not be complete. Incomplete Chain of Custody will result in delay of processing or inability to process as Chain of Custody. (required)

Custody Change <small>(required)</small>	Responsible Party <small>(required)</small>	Print Name (Last, First Middle) <small>(required)</small>	Sign <small>(required)</small>	Date (mm-dd-yyyy) <small>(required)</small>
Specimen provided for testing	Donor Full Name, or Legal Representative or Medical Staff (if donor unable to sign)			
Specimen received and sealed	Collector			

This section below will be filled out by Mayo Clinic Laboratories personnel upon receipt of the specimen. (required)

Custody Change <small>(required)</small>	Responsible Party <small>(required)</small>	Print Name (Last, First Middle) <small>(required)</small>	Sign <small>(required)</small>	Date (mm-dd-yyyy) <small>(required)</small>
Specimen received by Mayo Clinic Clinical and Forensic Toxicology Laboratory	Mayo Clinic Lab Personnel			



Ship specimens to:

Mayo Clinic Laboratories
3050 Superior Drive NW
Rochester, MN 55905

Customer Service: 800-533-1710

Billing Information

- An itemized invoice will be sent each month.
- Payment terms are net 30 days.

Call the Business Office with billing-related questions:
800-447-6424 (US and Canada)
507-266-5490 (outside the US)

Visit www.MayoClinicLabs.com for the most up-to-date test and shipping information, or scan the code to learn more.

To fill out this form electronically, scan the code. Then print and send the form with the specimen.

Patient Information (required)

Patient ID (Medical Record No.)	Patient Name (Last, First Middle)
Birth Date (mm-dd-yyyy)	Client Account No.

Mayo Clinic Laboratories does not perform workplace drug testing.

URINE TESTING

CDA7X Drug Abuse Survey with Confirmation, Panel 9, Chain of Custody, Random, Urine

PDSUX Drug Screen, Prescription/Over the Counter, Chain of Custody, Random, Urine

PANOX Pain Clinic Survey 10, Chain of Custody, Random, Urine

OXYSX Oxycodone Screen, Chain of Custody, Random, Urine

ADLTX Adulterants Survey, Chain of Custody, Random, Urine

INDIVIDUAL URINE DRUG TESTING

6MAMX 6-Monoacetylmorphine, Chain of Custody, Random, Urine

AMPHX Amphetamines Confirmation, Chain of Custody, Random, Urine

ETGX Ethyl Glucuronide Confirmation, Chain of Custody, Random, Urine

BARBX Barbiturates Confirmation, Chain of Custody, Random, Urine

BNZX Benzodiazepines Confirmation, Chain of Custody, Random, Urine

BUPMX Buprenorphine and Norbuprenorphine, Chain of Custody, Random, Urine

COKEX Cocaine and Metabolite Confirmation, Chain of Custody, Random, Urine

FENTX Fentanyl with Metabolite Confirmation, Chain of Custody, Random, Urine

MTDNX Methadone Confirmation, Chain of Custody, Random, Urine

OPATX Opiates Confirmation, Chain of Custody, Random, Urine

OXYCX Oxycodone with Metabolite Confirmation, Chain of Custody, Random, Urine

PCPX Phencyclidine Confirmation, Chain of Custody, Random, Urine

THCX Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC) Confirmation, Chain of Custody, Random, Urine

VLTUX Volatile Screen, Chain of Custody, Random, Urine

Urine temperature is within range of 90.5° F to 99.8° F
 Yes No Not measured
 If No, record temperature: _____

THCMX 11-nor-Delta-9-Tetrahydrocannabinol-9-Carboxylic Acid (Carboxy-THC) Confirmation, Chain of Custody, Meconium

COKMX Cocaine and Metabolite Confirmation, Chain of Custody, Meconium

OPTMX Opiate Confirmation, Chain of Custody, Meconium

PCPMX Phencyclidine (PCP) Confirmation, Chain of Custody, Meconium

ADDITIONAL TESTS (indicate Test ID and name)

BLOOD TESTING

VLTBX Volatile Screen, Chain of Custody, Blood

MECONIUM TESTING

DSM4X Drugs of Abuse Screen 4, Chain of Custody, Meconium

DSM5X Drugs of Abuse Screen 5, Chain of Custody, Meconium

MAMMX 6-Monoacetylmorphine (6-MAM) Confirmation, Chain of Custody, Meconium

AMPMX Amphetamine-Type Stimulants Confirmation, Chain of Custody, Meconium

REMARKS

Urine Collection Procedure

- Remove transport bag, collection cup, transport bottle, temperature strips and security tape from the Chain of Custody Kit.
- Complete the **Patient Information** and **Test** sections of this form.
- If this is an **unwitnessed collection**, do not allow the donor into restroom until steps below are followed. This will decrease the possibility for an adulterated collection.
 - Add bluing to the toilet water.
 - Tape the top of the toilet tank closed.
 - Secure sink and soap dispensers with tape.
 - Clear all areas, including cupboards and garbage containers, of agents that could be added to urine to void testing results.
 - Do not allow coats, bags, etc that could conceal adulteration material.
- Verify donor's identity. Instruct donor not to flush toilet or run water during collection. Give collection cup to donor and tell them to contribute a minimum of 60 mL of urine.

Following Collection

- Immediately verify sufficient volume was collected and check sample integrity by examining appearance, odor or other unusual characteristics.
- Immediately record temperature, if required.
- Have donor sign and date form under **Donor Signature/Date/Printed Name**. (If donor is not able to sign, a legal representative or medical staff member may provide their own Full Name/Sign/Date.)
- Pour a minimum of 30 mL** from the collection cup into the 60 mL urine transport bottle.

- Collection personnel sign or initial and date a security tape and seal over the cap of the specimen.
- Indicate 2 patient identifiers on the specimen container.** Either write these identifiers on the container or place a label with these identifiers on the container. Identifiers must match the information that was completed at the top of the form.
 - Patient name or unique patient identifier.
 - Patient number, hospital or specimen number, or birth date.
Do **not** place security tape over these identifiers.
- Indicate if a split specimen was collected or if there were collection problems.
- When specimens are sealed, print collector's name and sign and date form under **Collector Signature/Date/Printed Name**. **Note:** Date must match date on the security tape.
- Place specimen into rear pouch of transport bag and seal pouch by removing blue tape.
- Place **Chain of Custody form with original signatures** into the front pouch. Forward to Mayo Clinic Laboratories.
 - If blood, serum, or meconium is collected, complete steps 1, 2, 5, 7, 9, 10, 11, 12, 13, and 14.**
 - During multiple meconium collections, document each collection date, time, and name of collector in the remarks field.**
 - Blood or serum collection—Do NOT use alcohol or alcohol prep to clean arm before collection.**