



# Test Definition: FIL4S

Interleukin-4 (IL-4) Serum

## Overview

### Method Name

Multiplex array electrochemiluminescence

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

**Specimen Type:** Serum

**Collection Container/Tube:** Red or SST

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

### Specimen Minimum Volume

1 mL

### Reject Due To

Thawing:	Warm reject; Cold OK
----------	----------------------

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

## Clinical & Interpretive

### Reference Values

<0.5 pg/mL

**Performance****PDF Report**

No

**Day(s) Performed**

Monday, Wednesday, Friday

**Report Available**

5 to 11 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

[This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.](#)

**CPT Code Information**

83520

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FIL4S	IL-4, Serum	27161-9

Result ID	Test Result Name	Result LOINC® Value
FIL4S	IL-4, Serum	27161-9