



Test Definition: FABP2

Allergic Bronchopulmonary Aspergillosis Panel

II

Overview

Method Name

Enzyme Immunoassay (FEIA)
Ouchterlony Gel Immunodiffusion

NY State Available

No

Specimen

Specimen Type

Serum

Specimen Required

Collection Container/Tube:

Preferred: Red top

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL Serum

Collection Instructions:

1. Centrifuge and aliquot 1.5 mL of serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume

Serum: 1 mL

Reject Due To

Lipemia	Reject
---------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information

<https://www.eurofins-viracor.com/test-menu/>

Reference Values

Total IgE (IU/mL):

1-11 months: 0-12

1 year: 0-15

2 years: 1-29

3 years: 4-35

4 years: 2-33

5 years: 8-56

6 years: 3-95

7 years: 2-88

8 years: 5-71

9 years: 3-88

10 years: 7-110

11-14 years: 7-111

15-19 years: 6-96

20-30 years: 4-59

31-50 years: 5-79

51-80 years: 3-48

Aspergillus fumigatus IgE

<0.35 kU/L

Aspergillus fumigatus IgG

<46 mcg/mL

Aspergillus fumigatus Mix Gel Diffusion

Negative

Interpretation

Aspergillus fumigatus IgE

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal/Borderline
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	High Positive
4	17.50-49.99	Very High Positive
5	50.00-99.99	Very High Positive
6	>99.99	Very High Positive

Performance**PDF Report**

No

Day(s) Performed

Monday through Friday

Report Available

5 to 9 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86331

86001

86003

82785

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FABP2	Allergic Bronchopulmonary Panel II	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z3248	Total IgE	In Process
Z3249	Aspergillus fumigatus IgE	6025-1
Z3250	Class	102153-4
Z3251	Aspergillus fumigatus IgG	26954-8
Z3252	A. fumigatus Mix Gel Diffusion	23821-2