



# Test Definition: FMACR

Macroamylase

## Overview

### Method Name

Polyethylene Glycol Precipitation, Spectrophotometry

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

1 mL

### Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	NA
Icterus	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	30 days	

## Clinical & Interpretive

### Reference Values

Not Detected

**Performance****PDF Report**

No

**Day(s) Performed**

Tuesday, Thursday

**Report Available**

2 to 9 days

**Performing Laboratory Location**

Quest Diagnostics Nichols Institute

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**CPT Code Information**

82150

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FMACR	Macroamylase	15358-5

Result ID	Test Result Name	Result LOINC® Value
FMACR	Macroamylase	15358-5