



# Test Definition: FANTI

Anti-IgE

## Overview

### Method Name

ELISA

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

#### Collection Container/Tube:

**Preferred:** Red top

**Acceptable:** Serum gel

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL Serum

#### Collection Instructions:

1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerate

### Specimen Minimum Volume

Serum: 0.5 mL

### Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Ambient (preferred)	7 days	
	Refrigerated	7 days	
	Frozen	28 days	

## Clinical & Interpretive

### Clinical Information

---

Refer to [www.eurofins-viracor.com/test-menu](http://www.eurofins-viracor.com/test-menu)

**Reference Values**

Normal

**Performance****PDF Report**

No

**Day(s) Performed**

Wednesday

**Report Available**

7 to 11 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

83516

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FANTI	Anti-IgE	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FANTI	Anti-IgE	Not Provided