



# Test Definition: XSRM

Additional Sample for Reflex Oligoclonal Banding, Serum

## Overview

### Method Name

Only orderable as part of a profile. For more information see MSP3 / Multiple Sclerosis (MS) Cascade, Serum and Spinal Fluid.

Collection Only

### NY State Available

No

## Specimen

### Specimen Type

Serum

### Specimen Required

Only orderable as part of a profile. For more information see MSP3 / Multiple Sclerosis (MS) Cascade, Serum and Spinal Fluid.

### Collection Container/Tube:

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

### Collection Instructions:

1. Within 2 hours of collection, centrifuge and aliquot serum in plastic vial.
2. Label specimen as serum.

### Specimen Minimum Volume

0.5 mL

### Reject Due To

Gross hemolysis	OK
Gross lipemia	OK
Gross icterus	OK

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Ambient	14 days	
	Refrigerated	14 days	

## Clinical & Interpretive

## Performance

### PDF Report

No

### Specimen Retention Time

14 days

### Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Superior Drive

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

Not Applicable

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
XSRM	Additional sample for reflex OLIGS	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
XSRM	Additional sample for Reflex OLIGS	No LOINC Needed