



# Test Definition: FGLIP

Glipizide (Glucotrol)

## Overview

### Method Name

High Performance Liquid Chromatography with Ultraviolet Detection (HPLC-UV)

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Specimen Required

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 2 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

0.3 mL

### Reject Due To

|           |    |
|-----------|----|
| Hemolysis | NA |
| Lipemia   | NA |
| Icterus   | NA |
| Other     | NA |

### Specimen Stability Information

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Varies        | Refrigerated (preferred) | 14 days  |                   |
|               | Ambient                  | 72 hours |                   |
|               | Frozen                   | 180 days |                   |

**Clinical & Interpretive****Reference Values**

Units: ng/mL

Plasma insulin concentrations have been shown to increase only when plasma glipizide concentrations exceeded 200 ng/mL.

Toxic range has not been established.

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

7 to 10 days

**Specimen Retention Time**

2 weeks

**Performing Laboratory Location**

Medtox Laboratories, Inc.

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**CPT Code Information**

80299

**LOINC® Information**

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
|---------|-----------------|--------------------|

|       |                       |         |
|-------|-----------------------|---------|
| FGLIP | Glipizide (Glucotrol) | 10539-5 |
|-------|-----------------------|---------|

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| Z1180     | Glipizide        | 10539-5             |